

Vitrex[®] Flex 3 Depth Adjustable Safety lancet

Vitrex Flex 3 safety lancet is an All-in-One lancet for emergency use due to its adjustable design, which includes three different penetration depths for a single unit. The lancet is intended for blood glucose and other in vitro diagnostic tests.

USER CONVENIENCE: Pre-loaded and pressure activated design. The trigger position indicates the before and after use.

USER SAFETY: Needle is fully shielded before and after use.

COMFORT: High quality tri-bevel needle with high speed penetration.

PATIENT SAFETY: Single-use only. Sterilised by gamma irridiation.



Improved safety cap removal – cap automatically semi ejects, when turned to remove.



Avoid misfiring when adjusting the depth. Ergonomic design for improved holding.

Technical Specifications

Colour Code		Specification	Penetration Depth	Article No.
Blue	\bigcirc	23G	1.30 mm - 1.80 mm - 2.30 mm	395813
Yellow	\bigcirc	26G	1.30 mm - 1.80 mm - 2.30 mm	396813
Orange		28G	1.30 mm - 1.80 mm - 2.30 mm	397813

Shelf life 5 years



Packaging: Primary packaging: Transit packaging: Pallet:



Penetration depth adjustment is done without touching the activation button.